

Invacare Matrx Flo-tech Image

The *Invacare Matrx Flo-tech Image* cushion uses a combination of a slim, shaped foam cushion base and a large gel sac covering the sacral area, increasing even weight distribution and immersion of the ischial tuberosities. This design delivers a high degree of protection to users who are at High Risk of developing pressure ulcers. This slim cushion format is ideal for both chair and wheelchair use.

Slimline contoured foam and gel combination cushion

Some users prefer a slimline cushion, finding it more compatible with their daily living activities. The *Invacare Matrx Flo-tech Image* combines the benefit of a contoured foam upper surface with a large gel sac covering the ischial and sacral area, offering increased support and reducing the occurrence of peak pressures.

Contoured seating surface with gel sac

Offers improved stability and postural support; moulded foam castellations and

gel combine effectively to reduce shear and friction forces on tissue.

Features and Options



Ingress Resistant Cover

Cover is two-way stretch, waterresistant and vapour-permeable. Cover features an anti-slip base with carry handle.



Available in 23 sizes

The *Invacare Matrx Flo-tech Image* is available in a range of 23 sizes, from 36 x 33 to 46 x 41 cm (360 x 330 to 460 x 410 mm).



Integral Sag Compensator option

Option of sag compensator to accommodate the sag in wheelchair seats, and ensure the cushion is correctly presented.



Ultra-Fresh™ added to all *Flotech* foam for antimicrobial and odour protection

Antimicrobial treatment that has been added to all **Flo-tech** foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.