



Invacare **Matrix Flo-tech Contour Lo-Back Cushion**

The new *Invacare Matrix Flo-tech Contour Lo-Back Cushion* has been successfully designed to have a lower rear edge which will consequently allow the user to sit further back in the seat, effectively optimizing the support surface contact area for pressure distribution, postural support and reducing sacral area interface pressures.

Contour Lo-Back Cushion optimizes support surface contact area

Comparison with standard wheelchair cushion, showing improved posture through use of the **Matrix Flo-tech Contour Lo-Back Cushion**.

Moulded Foam Cushion

Moulding allows **Matrix Flo-Tech** cushions to provide a three-dimensional form, which cannot be emulated in cut

foam. Moulded foam is highly durable, increasing the product's longevity.

Features and Options



Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.



Available in 83 sizes

The *Invacare Matrix Flo-tech Contour Lo-Back Cushion* is available in a range of 83 sizes, from 36 x 33 to 60 x 55 cm (360 x 330 to 600 x 550 mm).



Contoured seating

Offers improved stability and postural support. Castellations reduce sheer and friction forces on tissue



Ultra-Fresh™ added to all Flo-tech foam for antimicrobial and odour protection

Antimicrobial treatment that has been added to all **Flo-tech** foam to help make it resistant to the growth of odour causing bacteria. Helps provide longer lasting freshness.